

Kit	Product Desc.	Part Number		Instrument																
				Stratagene Mx			ABI									Qiagen	Eppendorf Master cycler	BioRad		Roche
QPCR		400 rxns	4000 rxns	3000P	3005P	4000	7000	7300	7500	7500 Fast	7700	7900HT	Step One	Step One Plus	Rotor-Gene	ep realplex	CFX96	iCycler iQ	Light Cycler 480	
Probe																				
QPCR Master Mix	Brilliant II QPCR Master Mix	600804	600815	+++	+++	+++										++	+	+	+	++
	Brilliant II FAST QPCR Master Mix	600845	600846	+++	+++	+++										++	+	+	+	++
	Brilliant II QPCR Low ROX Master Mix	600806	600817	+++	+++	+++			+++											
	Brilliant II QPCR High ROX Master Mix	600805	600816				+++	+++			+++	++	+++							
	Brilliant III Ultra-Fast QPCR Master Mix	600880	600881	?	?	?				+++		+++		+++	+++	+++	++	+++		+++
One-Step QRT-PCR	Brilliant II QRT-PCR Master Mix, 1-step	600809	600818	+++	+++	+++										++	+		+	++
	Brilliant II QPCR Low ROX Master Mix, 1-step	600837	600841	+++	+++	+++			+++											
	Brilliant II QPCR High ROX Master Mix, 1-step	600838	600842				+++	+++			+++	++	+++							
	Brilliant III Ultra-Fast QRT-PCR Master Mix	600884	600885	?	?	?				+++		+++		+++	+++	++	+++		+++	
Two-Step QRT-PCR	Brilliant II QRT-PCR AffinityScript Master Mix, 2-step	600827		+++	+++	+++	++	++	++		++			++	++	+		+	++	
SYBR																				
QPCR Master Mix	Brilliant II SYBR Green QPCR Master Mix	600828	600831	+++	+++	+++										++	+	+	+	++
	Brilliant II FAST SYBR QPCR Master Mix	600843	600844	+++	+++	+++										++	+	+	+	++
	Brilliant II SYBR QPCR Low ROX Master Mix	600830	600833	+++	+++	+++			+++											
	Brilliant II SYBR QPCR High ROX Master Mix	600829	600832				+++	+++			+++	++	+++							
	Brilliant III Ultra-Fast SYBR Green QPCR Master Mix	600882	600883	?	?	?				+++		+++		+++	+++	++	+++		+++	
One-Step QRT-PCR	Brilliant II SYBR Green QRT-PCR Master Mix, 1-Step	600825	600826	+++	+++	+++										++	+		+	++
	Brilliant II QRT-PCR SYBR Low ROX Master Mix, 1-step	600835	600839	+++	+++	+++			+++											
	Brilliant II QRT-PCR SYBR High ROX Master Mix, 1-Step	600836	600840				+++	+++			+++	++	+++							
	Brilliant III Ultra-Fast SYBR Green QRT-PCR Master Mix	600886	600887	?	?	?				+++		+++		+++	+++	++	+++		+++	
Two-Step QRT-PCR	Brilliant II SYBR Green QRT-PCR AffinityScript Master Mix, 2-step	600834		+++	+++	+++	++	++	++		++			++	++	+		+	++	

+++	Best choice
++	Good choice
+	Fair choice
?	Pending validation